

The Richie Brace Application



Step 1

Open front and back straps of the Richie Brace.



Step 2

Slide brace into the shoe.



Step 3

Slip foot into shoe through the back side of the brace.



Step 4

Kick heel back into the shoe.



Step 5

Position leg straight over the top of foot and tight front straps.



Step 6

Feed the back strap through outside slot. Twist upright of brace for better angle.



Step 7

Back strap should be tightened as tight as possible.



Step 8

Tie shoe laces.

Distributed By



635 Executive Dr. | Troy MI, 48083
P 800-521-2192 | 248-588-7480
BeckerOrthopedic.com | f t in