


Model 1013

Free Motion Instructions:

To construct this joint as free motion, remove the lever screw and bushing and replace with the included ratcheting hole plugs. Insert a plug into both sides of the knee joint and apply pressure.



©2021 Becker Orthopedic Appliance Co.
All rights reserved.
Revision 05/01/21

 **BECKER**
635 Executive Dr. | Troy, MI 48083
P 800-521-2192 | 248-588-7480
BeckerOrthopedic.com | f t in

EC REP

Acorn Regulatory Consultancy
Services Limited
Knockmorris Cahir Co. Tipperary
Ireland, Postcode: E21 R766
P 012 4626 8456
F 012 4626 8648

