

Bail Lock Integrated Strap System (BLISS) ASSEMBLY INSTRUCTIONS

Model MX-003-BLISS

- Secure the proximal strap with the enclosed screws to the high section of the KAFO (B1). Loop the bail straps over each lever and insert the supplied plastic caps into the lever openings (B4).
- Adjust strap length as required (B3). Be sure to set the length of the adjustment strap to ensure locked position is achieved at full extension.



©2021 Becker Orthopedic Appliance Co.
All rights reserved.
Revision 05/01/21



635 Executive Dr | Troy, MI 48083
P: 800-521-2192 | 248-588-7480
BeckerOrthopedic.com | f t in

EC REP

Acorn Regulatory Consultancy
Services Limited
Knockmorris Cahir Co. Tipperary
Ireland, Postcode: E21 R766
P 012 4626 8456
F 012 4626 8648



Bail Lock Integrated Strap System (BLISS) ASSEMBLY INSTRUCTIONS

Model MX-003-BLISS

- Secure the proximal strap with the enclosed screws to the high section of the KAFO (B1). Loop the bail straps over each lever and insert the supplied plastic caps into the lever openings (B4).
- Adjust strap length as required (B3). Be sure to set the length of the adjustment strap to ensure locked position is achieved at full extension.



©2021 Becker Orthopedic Appliance Co.
All rights reserved.
Revision 05/01/21



635 Executive Dr | Troy, MI 48083
P: 800-521-2192 | 248-588-7480
BeckerOrthopedic.com | f t in

EC REP

Acorn Regulatory Consultancy
Services Limited
Knockmorris Cahir Co. Tipperary
Ireland, Postcode: E21 R766
P 012 4626 8456
F 012 4626 8648

